

Summary of Analytic Considerations and Test System Monitoring

IRMA TRUPOINT® Blood Analysis System	
Potential Cause of Analytic Concern	Automatic Control/Monitoring
Environmental Considerations	
• Sample/environment/cartridge is below temperature during use	Yes
• Sample/environment/cartridge is above temperature during use	Yes
• Too much ambient light during use	No Effect on Instrument
• Too little ambient light during use	No Effect on Instrument
• Too much humidity during use	Yes
• Too little humidity during use	Yes
• Too much atmospheric pressure during use	Yes
• Too little atmospheric pressure during use	Yes
• Too much airflow during use	Yes
• Too little airflow during use	Yes
• Too much vibration during use	Yes
• Dust laden environment during use	Yes
• Exposure to RFI/EMI during use	Yes
• Shipping and storage issues	Yes
Operator Considerations	
• Not enough specimen when running test	Per Lab QC Competency Program
• Too much specimen when running test	Yes
• No specimen when running test	Yes
• Specimen contaminated microbially	Yes
• Specimen contaminated chemically	Yes
• Specimen contaminated with tissue fluids	Per Lab QC Competency Program
• Specimen begins to clot prior to testing/non-heparinized sample used	Per Lab QC Competency Program
• Improper anticoagulation material used	Yes
• Syringe used to inject sample is too small	Yes
• Injection of sample too fast or too slow	Per Lab QC Competency Program
• Specimen has air bubbles	Per Lab QC Competency Program
• Wrong assay cartridge used for test	Yes
• Used assay cartridge used for test	Yes
• Specimen added to cartridge prematurely	Yes
• Specimen not added to cartridge in a timely fashion	Yes
• Assay cartridge not fully inserted into instrument	Yes
• QC is not performed at required interval(s)	Automatic with Lock-out Enabled
• Patient test performed without appropriate QC	Automatic with Lock-out Enabled
• Patient ID is not entered into instrument	Automatic with Lock-out Enabled
• Operator ID is not entered into instrument	Automatic with Lock-out Enabled
• Untrained operator uses instrument	Automatic with Lock-out Enabled
• Instrument used during in-transit situations	Yes
• Incorrect entry of lot code or calibration code into instrument	Yes
• Lot code or calibration code is not entered into instrument	Per Lab QC Competency Program
• Expired cartridge used	Yes
Analysis Considerations	
• Power not available to start testing	Yes
• Power not available to complete testing	Yes
• Proper assay is not recognized by instrument	Yes
• Instrument software malfunctions	Yes
• Sensors do not function properly during test	Yes
• Detector does not function during test	Yes
• No reagent or calibrant gel during test	Yes
• Test result is inconsistent with patient's condition	Per Lab QC Competency Program



LifeHealth, LLC
 2656 Patton Road, Roseville, MN 55113 USA
 T: 1-855-762-8378, *1-651-638-1000
www.lifehealthmed.com

This information is the property of LifeHealth, LLC. It is based on LifeHealth's internal evaluations, current interpretation of the CMS CLIA Interpretative Guidelines in effect at the time of its issuance, and information in the public domain. Its sole intended use is to provide LifeHealth's internal evaluation summary of products related to the CLIA Interpretative Guideline for Equivalent Quality Control Options. Any other use or conclusions derived from this information is disclaimed by LifeHealth.